

Versatile and Stylish



Versatility knows no bounds with the Colt Plus. Enjoy simple disassembly with one hand and a bright LED headlight package. You won't miss a beat with this scooter, which includes a convenient front storage basket. Anything is possible with the Colt Plus!



Specifications

Weight Capacity	25 st (158.75 kg)
Maximum Speed ¹	Up to 4 mph
Ground Clearance ³	57 mm (2.25")
Turning Radius ³	1370 mm (54")
Overall Length ³	1190 mm (47")
Overall Width ³	560 mm (22")
Max. Obstacle Climb	78.7 mm (3.1")
Drive Tyres	254 mm front & rear - pneumatic
Suspension	Limited
Total Weight Without Batteries ^{3, 4}	62.5 kg (137.7 lbs)
Heaviest Piece When Disassembled ³	Main Frame: 29.7 kg (65.4 lbs)
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical
Battery Requirements ³	2 x 32 Ah deep cycle
Battery Weight ⁴	-
Battery Charger	4-amp off-board
Colour	Candy Apple Red

Features

- Comfortable fold flat seat with armrests
- High visibility front and rear LED lighting
- Wrap around delta tiller
- Feather touch disassembly
- High level tiller charger
- Exclusive & stylish low profile tyres
- Convenient storage in tiller

Available in:

- Candy Apple Red



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.